OIL-FREE

AIR FOR LIFE

DK50 D
DK50 DM

INSTALLATION, OPERATION AND MAINTENANCE MANUAL

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CE 0434

CONTENTS

1. GENERAL INFORMATION.....	2
INDICATIONS FOR USE	2
OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY	2
MARKING.....	2
WARNINGS.....	2
General safety warnings	3
Safety warnings regarding the protection against electric current.....	3
ALERT NOTICES AND SYMBOLS	4
USAGE ON PURPOSE.....	5
STORAGE AND TRANSPORT CONDITIONS.....	5
2. DEVICE DESCRIPTION.....	6
3. TECHNICAL DATA.....	8
4. OPERATION	9
INSTALLATION AND FIRST PUTTING INTO OPERATION.....	9
Unfixing of the air pump	9
Compressed air connection	10
Electric connection	10
First putting into operation.....	11
Additive accessory	11
PERSONNEL	12
Switch on compressor.....	12
Running on compressor	12
Cleaning of air intake filter.....	13
Cooling failure alarm	13
Cleaning the compressor	13
5. MAINTENANCE	14
INFORMATION ON REPAIR SERVICE.....	14
INTERVALS OF MAINTENANCE	14
Check of safety valve	15
Replacement of filters of air pump.....	15
Replacement of filter in blow filter	16
Replacement of filter in pressure regulator	16
Checking of capacity	17
Fixing of compressor before shipping	17
CUT-OFF.....	17
DISPOSAL OF THE APPLIANCE	17
6. TROUBLESHOOTING	18
7. SPARE PARTS.....	19
8. ELECTRIC AND PNEUMATIC DIAGRAMS	20
WIRING DIAGRAM	20
PNEUMATIC DIAGRAM.....	21
9. PARTS LIST	64
10.WARRANTY	67

GENERAL INFORMATION

1. GENERAL INFORMATION

INDICATIONS FOR USE

The EKOM DK50 D, DK50 DM is indicated for use as medical air compressor to supply a source of clean, oil-free pressurized air for use with medical ventilators.

OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY

The installation, operation and maintenance manual is an integral part of the appliance. It is necessary to always keep this document close to the appliance. Strict observance of this manual is a prerequisite for the correct operation of the appliance.

Rx only

USA Federal law restrict this device to sale by or on the order of a physician.

MARKING

Products marked with CE mark of compliance meet safety requirements of European Union (93/42/EEC).

WARNINGS

- The safety of operating personnel and failure-free operation of the appliance are ensured only when using the original components of the appliance. It is possible to use accessories and spare parts mentioned in the technical documentation or clearly permitted by the manufacturer.
- When used with non authorized accessories or consumable material, the manufacturer cannot assume responsibility for the safe operation and functionality of the device.
- The warranty does not cover damages that originate due to the use of non authorized accessories or consumable material other than those recommended by the manufacturer.
- The manufacturer assumes responsibility regarding safety, reliability and function only if
 - the installation, calibration, amendments, extensions and repairs are made by the manufacturer or his representative or a service organization authorized by the manufacturer,
 - the appliance is used in accordance with the installation, operation and maintenance manual.
- The installation, operation and maintenance manual, at the time of printing, correspond to the design of the appliance and its state according to the relevant safety and technical standards. The manufacturer reserves all copyrights for the given wiring schemes, methods and names.

General safety warnings

The manufacturer developed and constructed the appliance so that damage would not occur when the appliance is used for its intended purpose. The manufacturer considers it his obligation to describe the following safety measures in order to avoid further damages.

- When operating the appliance, it is necessary to observe laws and regional regulations valid in the place of usage. In order to ensure safe course of works, the operator and user are responsible for the observation of regulations.
- The original packaging should be kept for the possible return of the unit. Only original packaging guarantees an optimal protection of the appliance during transportation. If it would be necessary to return the appliance during warranty period, the manufacturer is not responsible for damages caused by incorrect packing.
- If an undesirable event occurs in the operation of appliance, the user is obliged to immediately inform his supplier to this event.
- This product is not intended for use in areas with the risk of explosion. Possible explosion hazard if the compressor is used in the presence of flammable anesthetics.
- Never put oxygen or nitrous oxide into the compressor. The electrical components are not approved for oxygen or nitrous oxide use.
- This device can be used only with ventilator equipped with low pressure alarm.

Safety warnings regarding the protection against electric current

- The appliance must only be connected to an appropriate power source that has correct grounding. Reliable earth ground can only be achieved when equipment is connected to a receptacle marked "hospital grade".
- Prior the connecting the compressor, verify whether the mains voltage and frequency specified on the apparatus are in accordance with the local supply.
- Prior to putting into operation, check for possible damages on the appliance and the air connectors. Damaged cables and sockets/plugs must be replaced immediately.
- In the case of a dangerous situation or a technical failure, immediately disconnect the appliance from mains supply.
- During all repairs and maintenance:
 - ensure that the mains plug is removed from the power socket,
 - pressure pipes must be air vented,
 - pressure must be released from pressure tank.
- Only a qualified expert can install this appliance.

GENERAL INFORMATION

ALERT NOTICES AND SYMBOLS

For your information, the symbols below are used in the installation, operation and maintenance manual, on packaging materials and on the product:

Attention, see instructions for use

Caution, risk of electric shock

Consult instructions for use

CE mark of compliance

Caution, hot surface

Handling mark on packaging – Fragile, handle with care

Handling mark on packaging – This way up (Vertical position of freight)

Handling mark on packaging – Protect against humidity

Mark on packaging – Recyclable material

Handling mark on packaging – Temperature of storage and transport

Handling mark on packaging – Limited stacking

Alternating current

Protective Earth (Ground)

Equipotentiality

Fuse

USAGE ON PURPOSE

- The appliance can be installed and operated only in dry, ventilated and dust-free area. Climatic conditions for operation - see Technical data.
- The compressor must stand on a flat, sufficiently stable base with.
- The compressor must not be exposed to rain. The appliance must not be used in humid or wet environments. Furthermore, it is forbidden to use it in the presence of flammable liquids or gasses.
- Prior to the connection of compressor to the respiratory appliance, it is necessary to assess whether the available medium meets the requirements of the given appliance of use. Please, refer to the Technical data for this purpose.
- Usage, other than the compressors intended use is not considered to be safe. The manufacturer is not responsible for the damages implied if the compressor is used for other uses. Such a risk is assumed only by operator/user.

STORAGE AND TRANSPORT CONDITIONS

The compressor is dispatched from the factory in transport packing with the pump fixed. Thus the appliance is protected from damage during transport. P

During transport, if it is possible, always use the original packaging of the compressor. Always transport the compressor up right.

During transport and storage, protect the compressor against humidity, contamination and extreme temperatures. Compressors that have original packaging can be stored in warm, dry and dust-free areas.

Keep the packaging material, if possible. If it is not possible, dispose of the packaging material in an environment-friendly way.

The compressor can only be transported pressure-free. It is necessary to release the pressure from the pressure tank and pressure hoses prior to transport, and to release any possible condensate. Secure motor inside before shipping!

Prior to transport it is necessary to fix motor inside the compressor (Chapter 5.)

DEVICE DESCRIPTION

2. DEVICE DESCRIPTION


1. PRESSURE - the output air pressure manometer (with accuracy $\pm 5\%$)
2. HUMIDITY - the humidity indicator (performance)
3. Cooling failure indicator
4. Power indicator
5. Hour meter
6. Main input, power switch, primary fuses
7. Equipotential pin
8. Socket for the supply lead
9. WALL - input air coming from outside pressure air source
10. OUT - output air from compressor
11. Pressure tank
12. Air intake filter
13. Pressure regulator with filter
14. Filter and water trap
15. Cooler
16. Control electronics
17. Safety valve
18. Internal back cover
19. Tube for supporting equipment

The appliance consists of an oil-free piston compressor driven by a single-phase electric motor, of filtering and drying equipment, of a pressure tank and of a pressure regulator. The compressor is placed in a box lined with noise reducing material. The appliance thus produces filtered compressed air without any traces of oil with reduced humidity contents.

The air compressor has an input, air fitting labeled WALL (9) that allows you to connect to an outside air source. When an outside compressed air source that is connected to the inlet WALL fitting supplies air above adjusted output pressure the air pump will not be activated even though the power switch is in the ON position. Instead, the air from the outside source will be directed to the compressor's air output. When the outside air source pressure is less than adjusted output pressure (or is not connected) and the power switch is in the ON position, the air pump is activated and supplies compressed air through the air output.

The air pump sucks the atmospheric air through air intake filter (12). Compressed air is cooled in the cooler (15), continue in the filter (14) and condensed liquid from the water trap is automatically released to the evaporative tray. The air continues through the check valve to the pressure tank (11) and through the pressure regulator with the output filter (13) to the OUT (10). The fans operate while the air pump works, or also during the pause, when the temperature in the cabinet is over 40°C (104°F). Any increase of temperature in the cabinet over 80°C (176°F) - cooling failure - is signalled by the red indicator (3) and by acoustic alarm. The Control electronics (16) control running the air pump according to pressure in the pressure tank (11). A safety valve (17) prevents the maximal allowable pressure in the pressure tank.

DEVICE DESCRIPTION

TECHNICAL DATA

3. TECHNICAL DATA

TYPE	DK50 D				DK50 DM								
Output flow at pressure 3.5 bar (51 psig)	Lit.min ⁻¹				40				60				
Output flow max. at pressure 3.5 bar (51 psig)	200 Lit.min ⁻¹ for 2 s *												
Rated voltage / Frequency	V/Hz	230/50	230/60	USA Model 120/60	100 /50-60	230/50	230/60	USA Model 120/60					
Nominal current	A	2.8	2.8	5.6	5,6	3.4	4.3	8.6					
Filtration of air to	µm	5											
Dew point depression at 40 Lit.min ⁻¹ , 3 bar (43.5 psig), 20°C (68°F)	5°C (9°F) under the ambient temperature												
Outlet connection	RECTUS 21KS KS08 MPC At pleasure DISS 1160-A (3/4"-16 UNF)												
Sound level	dB(A)	≤ 49	≤ 51	≤ 51	≤ 51	≤ 50	≤ 52	≤ 52					
Mode of the operation	Continuous - S1												
Release of condensed water	Automatic												
Indication of cooling failure	Acoustic and optical by 80°C (176°F)												
Output pressure	3.0 bar (43 psi) Adjustable to max. 3.5 bar (51 psi)												
Cut-in pressure (backup)	3.0 bar (43 psig) Adjustable to max. 3.5 bar (51 psi)												
Pressure tank capacity	Lit	5											
Working pressure of the pump	5.0 - 7.0 bar (72.5 - 101.5 psi)												
Operation pressure of safety valve	8.0 bar (116 psi)												
Adjustment of the output air pressure	Pressure regulator												
Dimensions of compressor w x d x h	400x430x780 mm (16x17x31 in)				400x500x780 mm (16x20x31 in)								
Dimensions of packaging w x d x h	480x530x835 mm (19x20x33 in)				480x580x835 mm (19x23x33 in)								
Net weight	46 kg (101 lbs)				59 kg (130 lbs)								
Gross weight	60 kg (132 lbs)				73 kg (161 lbs)								
Implementation according EN 60601-1, EN 12021	Appliance of type B class I.												
Classification acc. to MDD 93/42 EEC, 2007/47 EC	II b												

* With 0.6 bar (8.7 psi) drop

Climatic conditions for storage and transport

Temperature -25°C to +55°C (-13°F to +131°F), 24 hrs +70°C (+158°F)

Relative air humidity 10% to 90 % (no condensation)

Relative air humidity for seaworthy packing 10% to 100% (included condensation)

Climatic conditions for operation

Temperature +5°C to +40°C (+41°F to +104°F)

Relative air humidity up to +95%

IPX0 Rating

4. OPERATION

INSTALLATION AND FIRST PUTTING INTO OPERATION

It's forbidden to use compressor after unpacking as far as it doesn't adjust ambient temperature.

The compressor must be installed and put into operation for the first time by a qualified expert. His obligation is to also train the operating personnel about the usage and common maintenance of the appliance. Installation and training of personnel should be acknowledged by his signature in the certificate of installation.

Prior to first putting into operation, all securing fixings serving for the fixation of the appliance during transport must be removed. In the case of switching the compressor on without the prior removal of fixings, there will be a risk of permanent damage to the product!

The compressor does not contain an alternate power supply (backup supply).

This device can be used only with ventilator equipped with low pressure alarm.

It is forbidden obstruct the air intake filter and the ventilating outlets on the lower part of appliance. Forbidden to lift unit by handle! Handle is appointed for moving unbraked product on the floor.

If the compressor is equipped with main source of air, standby source of air must be available.

Unfixing of the air pump

After unpacking, place the compressor upright with the castor wheels firmly on the ground.

- Remove the front and back cover unscrewing screws 6 pcs on the side of cabinet.

OPERATION

- Unscrew the four screws and remove the internal back cover (18).
- Remove all fixing elements (21) above and under the air pump (20).
- Screw the internal back cover, connect the ground wire, screw the front and back cover and insert the air intake filter.

Please retain the fixing elements for future transportation of the compressor.

Compressed air connection

The medical compressor is equipped with fast-on couplings **WALL** (9) and **OUT** (10) to the rear of the cabinet.

To the fast-on coupling **OUT** (10) - output compressed air - connect the pressure hose leading to the applicable appliance / respiratory apparatus.

To the fast-on coupling **WALL** (9) - input compressed air - connect pressurized air from a central distribution system, if necessary. The air of this distribution system is automatically connected through the compressor to the outlet OUT port. In such a connection, the compressor serves as back up source of the compressed air. In the case of the reduction of the pressure of air from central distribution system under adjusted output pressure, compressor is automatically switched on and there is continuously pressurized air at the output from compressor.

Please note that the air supply from the central distribution entering the compressor must be medical grade air (dirt size, humidity). The compressor not modify the air from central distribution!

The hose connecting the compressor to the breathing apparatus must not pass through a cold environment i.e. placed on the ground. It should be as short as possible with no kinks (it may have impact on condensation of water inside the hose).

Electric connection

The compressor is delivered with a plug containing an appropriate protective contact. It is essential to respect the local electrical regulations. The mains voltage and frequency must comply with the information on the Data Label.

The supply lead must not be stressed or have any tension exerted upon it, and must always be free.

- The socket must be accessible due to safety reasons, in order for the appliance to be safely disconnected from the mains supply in the case of danger.
- The relevant current circuit must be safeguarded.
- Connection of the equipotential connection (7) with other appliances must be completed according to local valid electrical regulations.
- Fasten down the supply lead through the socket (8).

First putting into operation

- Check whether all fixing screws used during transport were removed.
- Check that the connections of the compressed air supplies are correct.
- Check for proper connection to the mains supply.
- Switch on the pressure switch (6) to position "I".

After first putting the compressor into operation, the pressure tank of compressor is pressurized to a 'switching-off' pressure and the compressor will then be automatically switched off. During the following time period the compressor works in an automatic mode, i.e. according to the consumption of compressed air, the compressor is switched on and off by the pressure switch.

During the operation, the appliance releases condensed water from the condenser (filter) automatically to the evaporative tray.

Additive accessory

Tube for supporting equipment - art. Nr. 604011346

Tube (19) is suitable for gripping of respiratory apparatus with accessories. Please remove the plug in the middle of the upper side before installation. Screw the tube into the socket and tighten it fast.

Technical data:

diameter 25 mm
height 1100 mm above the floor


Max. tube loading accepted is 18 kg. The centre of gravity of supporting equipment must be in the axis of the tube, max. 1100 mm above the floor. Other connecting consults with manufacturer. Manufacturer recommends supporting equipment disassemble before transport.

Set for higher relative air humidity - art. Nr. 604011497

1. Hose	PUR 6/4-130	062000337	1 pc
2. Hose	PUR 6/4-700	062000357	1 pc
3. Bottle	LUH-1	062000333	1 pc
4. Holder	3KB-972	023000564	1 pc
5. Washer	D5	043000002, 043000005	4 pcs
6. Screw	M5	041000131	4 pcs
7. Reduction	10/6	062000339	1 pc
8. Connector	L6/6	025400153, 025400154	1 pc
9. Clamp	self-adhesive	033400011	2 pcs

Dismount hose (A) and plug (B). Install reduction (7), connector (8) and plug-in hose (1). Screw holder (4), plug-in hose (2) and install bottle (3). Fix hose with two clamps (9).

OPERATION

Empty bottle regularly.

PERSONNEL

In the case of danger, switch the appliance off at the switch and pull out the mains plug.

Automatic operation. When the pressure in pressure tank falls down to a 'switch on' pressure, the compressor is automatically switched on. The compressor is automatically switched off when the pressure increases to the 'switching off' pressure in the pressure tank.

Switch on compressor

The compressor is switched on by mains switch (6) by putting it to position „I“. Switching on of the compressor is signaled by the green indicator (4).

Running on compressor

When the pressure of the air at the 'Inlet' port is reduced or if there is a reduction of pressure in the pressure tank, the compressor is switched on and the pressure tank is filled with compressed air. After the 'switching-off' pressure is reached, the compressor is automatically switched off. Pressure of the output air from compressor indicates the manometer (1).

During the operation of the compressor, the pointer of the humidity indicator (2) should be in the green zone. If the consumption of air from an appliance is larger than the maximal outflow, the pointer of indicator will move from the green zone and the air from compressor will not reach the declarable parameters (see Technical data).

The hour meter (5) indicates total running time of the air pump.

Cleaning of air intake filter

At least, once a week remove and clean the air intake filter (12) located on the front side. Wash the filter in warm soapy water, rinse thoroughly and allow it to dry. Insert the clean filter until the sucking hole for filter is not perfectly covered by filter.

Cooling failure alarm

Problems with the air pump cooling signalize the red indicator (3) and the acoustic alarm.

In the event of cooling failure alarm is necessary ensure alternate air pressure source and disconnect the compressor.

Cleaning the compressor

The product can be cleaned by the usual detergents. Please do not clean the product using a solution containing abrasive components, chemical solvents or another aggressive chemical.

MAINTENANCE

5. MAINTENANCE

INFORMATION ON REPAIR SERVICE

The guarantee and extended guarantee repairs are to be completed by the manufacturer, or an organization authorized by the manufacturer.

The manufacturer reserves the right to modify the appliance in a way that will not impact substantially on the operation of the appliance.

Only a qualified expert or the Customer Service Dept. of the manufacturer must perform repairs that exceed the framework of common maintenance. Use only spare parts and accessories approved by the manufacturer.

Prior to any maintenance or repairs, it is necessary to switch off the compressor and disconnect it from mains (pull out mains plug). Remove the front and back covers unscrewing screws on the side of the cabinet; disconnect ground connecting on the cover.

INTERVALS OF MAINTENANCE


Maintenance range	Page	Time interval	To be performed by
Cleaning and replacement of air intake filter	13	At least once a week, replacement after 4000 hours	personnel
Check of safety valve	15	1 x per annum	qualified expert
Replacement of filters of air pump (A, B)	15	2000 hours running	qualified expert
Replacement of filter in blow filter (C)	16	2000 hours running	qualified expert
Replacement of filter in the pressure regulator (D)	16	2000 hours running	qualified expert
Testing of the tightness of joints and inspection of appliance	Service documentation	2000 hours running or after two years	qualified expert
Replacement of piston complete with piston rings, O-rings	Service documentation	4000 hours	qualified expert
Bearing replacement	Service documentation	4000 hours	qualified expert

Check of safety valve

The safety valve is adjusted to 8 bars (116 psig) by the manufacturer, tested and marked. It must not be readjusted.

- Turn the screw of the safety valve (17) several rotations to left until safety valve puffs.
- Let safety valve freely puff only **for a short period of time**.
- Turn screw to the right until the limit is reached; the valve must be closed now.


The safety valve must not be used for depressurizing the pressure tank. It could endanger the function of the safety valve.

Replacement of filters of air pump

Caution – hot surface! About disassembly the internal back cover can be hot the internal surface.

Disassembly the internal back cover. These filter inserts are in the top cover of crank case (20).

- Using screwdriver, remove the flexible guard.
- Remove the used and contaminated filters (A), (B).
- Insert a new set of input filter (B).
- Place prefilter (A) into the flexible guard (by plain side into the crank case) and then to the opening.

MAINTENANCE

Replacement of filter in blow filter


- Loosen lock on filter (14) by pulling it down.
- Slightly turn the cover o filter and pull it out.
- Unscrew the holder of filter.
- Replace filter (C), screw the holder of filter.
- Place the cover of filter and secure it by turning to right till locking.

Replacement of filter in pressure regulator


Prior the intervention to the appliance it is necessary to decrease the pressure of air in pressure tank to zero and to disconnect the appliance from mains.

- Unscrew the cover of filter regulator (13) and pull it out.
- Unscrew the holder of filter.
- Replace filter (D), screw the holder of filter.
- Screw the cover of filter regulator.

Set up output pressure air

Switch on the compressor and set the output pressure air flow until the humidity indicator (2) remains in the green area. The control button of the filter regulator (13) unlock, rotate to set the demand output pressure. After setting the pressure, lock the control button of the regulator pushing down.

Checking of capacity

Before checking of capacity it is necessary to switch off compressor by power switch, let pressure air out of pressure tank and unlink output hose from fitting OUT (10). Next (Than) to switch compressor by power switch and measure (record) the time before points start - stop (2). The capacity of compressor is satisfied when the recorded time is not higher than 45 sec for DK50 D or 30 sec for DK50 DM.

Fixing of compressor before shipping

Before shipping it is required to fix the compressor to prevent movement. Procedure is illustrated on the page 9.

- Remove the front and back covers unscrewing screws on the side of the cabinet, disconnect the ground connector on the covers.
- Disassembly the internal back cover (18) unscrewing four screws.
- Insert fixing elements (21) under and up the air pump.
- Assembly the internal back cover (18), connect the ground connector on the covers, insert the covers and the air intake filter.

CUT-OFF

If the compressor is not going to be used for a long period of time, it is recommended that the appliance be disconnected from the mains supply and the air pressure in the pressure tank be released.

DISPOSAL OF THE APPLIANCE

- Disconnect the appliance from the mains supply.
- Release the air pressure in pressure tank.
- Dispose of the appliance according to local valid regulations.

The parts used in this product, when disposed, do not have negative impact on the environment.

TROUBLESHOOTING

6. TROUBLESHOOTING

Prior to repairing the appliance it is necessary to reduce the pressure of the air in the pressure tank to zero and disconnect the appliance from the mains supply.

Trained service personnel can only perform the activities connected to the troubleshooting guide.

PROBLEM	POTENTIAL CAUSE	TROUBLESHOOTING
Compressor does not start	No mains voltage	Switched off mains breaker in distribution system
	No current in pressure switch	Check Supply
		Replace the defective fuse
		Loose wire terminal – tighten
	Power cord defective – replace	
	Interrupted winding of motor, damaged thermal protection	Replace motor
	Defective capacitor	Replace capacitor
The humidity indicator is out of green zone	Jammed piston or another parts	Replace damaged parts
	Press. switch does not switch on	Change pressure switch
Compressor often switches on	Leakage of air from pneumatic distribution system	Control of pneumatic distribution system – seal loose connection
	Large consumption of compressed air	Do not exceed max. flow (see Technical data)
The operation of compressor is extended	Leakage of air from pneumatic distribution system	Control of pneumatic distribution system – seal loose connection
	Leakage of check valve	Clean check valve, replace sealing, replace check valve
	Worn out piston rings	Replace worn out piston rings
Red indication lights – state of overheating	Dirtily filter inserts	Replace the dirty inserts
	Dirty air intake filter	Clean or replace the dirty air intake filter
	Dirtily filter inserts	Replace the dirty inserts with new ones
	Cooling fans don't work	Replace fans or faulty thermo - switch
Compressor is noisy (knocking, metal noises)	Unit is hot, unventilated area	Relocate unit
	Damaged piston bearing or bearing on motor	Replace damaged bearing
	Loose (cracked) belt of the air pump hanger	Replace damaged hanger
Water comming out of outlets	Functionless water trap (13) or (14)	Clean or replace the water trap

7. SPARE PARTS


- Air intake filter (12) 4KB-552 025000011
- Set of filter inserts (A+B) 4EA-049 604041049
- Filter (C) 111511-40B 025200070
- Filter (D) AF20P-060S 025200113
- Fuse
- for DK50 D 230V T6,3A 038100004
- for DK50 D 120V, DK50 DM 230V
038100005 T10A
- for DK50 DM 120V T16A 038100007
- Insertion
- RECTUS 21SSTF08MPC 025500131
- DISS 1160-A 024000261

ELECTRIC AND PNEUMATIC DIAGRAMS


8. ELECTRIC AND PNEUMATIC DIAGRAMS

WIRING DIAGRAM

DK50 D

DK50 DM


Note

Wire entering the terminal should be marked with marking element with corresponding number.

Des.	Name	Des.	Name
HL1,HL2	Glowlamp	M	Electromotor
EV1,EV2	Ventilator	ST1,ST2	Thermo switch
PH	Hour meter	XC	Socket with fuses
YV	Solenoid valve	Cb1	Condensator

PNEUMATIC DIAGRAM

1. Intake filter
2. Input filter
3. Compressor
4. Solenoid valve
5. Sound absorber
6. Cooler
7. Microfilter with the water trap
8. Check valve
9. Pressure tank
10. Pressure regulator with filter
11. Air Outlet
12. Air Inlet
13. Main input, power switch, primary fuses
14. Pressure sensor
15. Humidity indicator
16. Manometer
17. Pressure relief valve
18. Fan
19. Water tray
20. Shuttle valve

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Extra equipment

This additional item is not sold to the subject of basic delivery, and must be ordered separately.

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Packing of basic equipment checked
by

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CERTIFICATE

The Product is manufactured according to production documentation and approved technical conditions. The product is complete and meets all prescribed conditions.

WARRANTY

The Manufacturer provides a warranty of the product, in duration of 24 months from the selling date, when the terms of warranty document are followed.

- In the case of a claim for repairs under warranty, a correctly completed guarantee document and certificate of quality and completeness of product must be submitted together with the product itself.
- During the warranty period all failures that occur to the product due to production faults will be repaired free of charge.
- The warranty does not apply to the product damaged during transportation or during incorrect storage.
- The warranty does not apply to failures due to incorrect operation of the device.
- The warranty does not apply to failures due to the use of the product for any other purposes other than which are described in the directions for use or as those which have been agreed with the manufacturer.
- The warranty does not apply to the product when any unauthorised modifications have been made.

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