MONZUN M1a

SUŠIČ S RIADENOU REGENERÁCIOU

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

Schéma zapojenia DK 50 2V / M1a

TABLE OF CONTENTS

1. INTENT AND USE	8
2. PRODUCT DESCRIPTION AND FUNCTION	8
2.1. Drying.....	8
2.2. Regeneration	8
3. RANGE OF DELIVERY	8
4. TECHNICAL DATA.....	8
5. ASSEMBLY AND INSTALLATION.....	8
6. PRODUCT OPERATION	8
7. PRODUCT MAINTENANCE	10
7.1. Replacement of filter insert	10
8. TROUBLESHOOTING	11
9. CONNECTION DIAGRAM.....	11
10. ELECTRIC WIRING DIAGRAMS	12
13. GUARANTEE.....	40

1. INTENT AND USE

The adsorption air drier **MONZUN** type **M1a** with a controlled regeneration is intended for drying and filtration of a pressure air of the oilfree compressors used for supplying the dental sets and devices in dental laboratories. The equipment referred to can also be used in other areas of medical practice, in laboratories, industry and places wherein a clean and dry air is necessary.

2. PRODUCT DESCRIPTION AND FUNCTION

The adsorption drier **MONZUN – M1a** works automatically in two phases :

2.1. Drying

Compressed air from the compressor (10) passes the cooler with ventilator (1) via adsorption chamber (7) (see diagram in chapter 8). By passing the chamber, the air is dried and filtered in the output filter (5 □m) situated in the upper part under a stopper (2) and is advancing further via the humidity sensor chamber (3) through the return valve (4) into the air chamber (6).

2.2. Regeneration

After each switching the compressor off by means of a pressure-operated switch, the solenoid valve (8) opens and will stay in the open status during the idle time of the compressor. Then a decompression of the adsorption chamber takes place, the adsorption substance is regenerated and accumulated water from its bottom part is expelled via a solenoid valve into the condensate vessel (9). When the humidity sensor evaluates the air passing the drier as „insufficiently dry“, it opens the solenoid valve of regeneration (5), through which it lets the regenerating air pass from the air chamber via the adsorption chamber within the time, unless the adsorption substance is regenerated to the necessary extent. The regeneration of adsorption substance only takes place during the time-out of the compressor operation.

The amount of air for the regeneration is different for the compressor DK 50 2V and for DK 50 PLUS.

For a long time reliable function of the drier, it is necessary to ensure its intermittent operation in a specified mode!

3. RANGE OF DELIVERY

Adsorption drier - type.....	M1a.....	1 pc
Instructions for use.....	NP-M1a-GB.....	1 pc
Filter intake 4KB-650.....	024000653.....	1 pc

4. TECHNICAL DATA

• rated voltage / frequency:	230 V ± 10% / 50 Hz
• rated power input:	35 VA
• amount of dried air:	140 lit.min ⁻¹ at 5 bar pressure
• mode of operation	intermittent- S3 - 60%
• working pressure:	5 - 7 bar
• drying stage:	atmospheric dew point -20°C
• filtration:	5 µm
• weight:	10 kg
• dimensions (height x width x depth):	700x345x200 mm

5. ASSEMBLY AND INSTALLATION

The equipment is factory-mounted to a corresponding compressor. The installation of a drier to compressors made so far should be consulted with the manufacturer.

6. PRODUCT OPERATION

A correct function of the drier depends on the compressor's operation and no attendance is required. The pressure vessel need not be sludged, because the pressure air entering the air chamber is already dried.

For a correct function of the drier it is necessary to proceed as follows:

- Observe an intermittent operation of the compressor in the 60% mode. The time of compressor operation should not, however, last longer than 10 minutes.
- Operate the compressor so that it works within the range of the on/off pressure of a pressure-operated switch set by the manufacturer. When the compressor works at a lower working pressure than that of

- the switching pressure, this is a sign of its overloading - high air offtake by an appliance, leakages in pneumatic distributions, aggregate or drier defects (see chapter Troubleshooting).
- Let the compressor connected to the electric network (neither switch the pressure-operated switch off, nor disconnect the power cord) – if the humidity sensor evaluates the air as „insufficiently dry“, it starts regenerating the drying chamber filling, that means due to the air passage from the air chamber, the pressure therein drops and the compressor can consecutively switch on several times. The compressor function will stop automatically after the drier regeneration.
 - In case that there is no air consumption and the regeneration takes place for more than 1 hour, make sure that there is an actual regeneration involved (check whether the air is withdrawn via the tubing into the condensate vessel) or the compressor or drier failures have occurred (see chapter Troubleshooting).


W a r n i n g !

After new installing of compressor with dryer M1a into operation (or after its longer inactivity, for example after a weekend, holidays, etc.) and its connection to electrical power network, the dryer regeneration starts. The compressor is put into operation automatically also without pressure air take-off by the appliance, and thus increase pressure in the air jet into cut-off pressure (see chapter Technical data of the Manual) when it turns off. Then the dryer regeneration takes place (passing of the air from air jet through dryer chamber), which can be expressed as a soft air hissing coming out of solenoid valve in the bottom part of the dryer. The air jet pressure decreases to closing pressure, the compressor turns on and increases pressure into air jet to cut-off pressure and turns off. The process of turning on and off of the compressor repeats in this way until the dryer is not regenerated sufficiently. This level is controlled in the dryer through the built-in humidity monitor – hygrostat.

The regeneration process can take several minutes (5 – 15 min) – for the new or previously already regenerated dryer during the previous operation, or several tens of minutes (30 – 120 min) – for the dryer being previously in its operation „overloaded“ with water vapour (for example due to the activity regime of the compressor outside the permitted extent, during the work in the environment with high relative humidity, etc.).

After finishing regeneration the whole process is stopped automatically.

Used jets: DK50 – 2V 0,7 mm
DK50 PLUS 0,5 mm

Note!

The mentioned values of the regeneration times are valid for compressor activity regime at closed outlet valve, which means without pressure air take-off by the appliance. In case of air take-off the mentioned times are prolonged.

In case the compressor does not finish the automatic dryer regeneration activity within 120 min., it is necessary to contact a supplier and request an opinion to operation.

For information the table shows the times of compressor regeneration cycles with dryer M1a on condition that air is not abstracted from compressor by appliance.

Compressor	Jet size	Closing pressure - cut-off pressure	Compressor activity time	Compressor break time - dryer regeneration
DK50 PLUS / M1a, DK50 PLUS S / M1a,	0.5 mm	4.5 – 6.5 bar	About 60 – 70 s	about 180 – 210 s
DK50 2V / M1a, DK50 2V S / M1a	0.7 mm	5 – 7 bar	about 30 – 40 s	About 90 – 110 s

In case of deviation from described activity please verify:

- Connection of the solenoid valve – compare arrow orientation on the valve body, which also marks correct air flow direction when regenerating
- Possible leakage via compressor untightness – close the outlet valve on the compressor, turn on the compressor and leave it in operation until it turns off at cut-off pressure. Turn off the circuit breaker and watch pressure on manometer, which means the pressure in air jet. The pressure decrease can not be bigger than 0.2 bar in 2 hours.
- The use of the correct jet – a jet is mounted between solenoid valve and opening (see Table). The dimension of an opening is marked on the jet body – number 7 equals 0.7 mm, number 5 equals 0.5 mm. When assembling it is important to keep jet orientation – aim the surface with a screwed cone at the valve.


7. PRODUCT MAINTENANCE

Warning!

Before making any interventions into the equipment, it is necessary to reduce the air pressure in the air chamber to zero and to disconnect the equipment from the electric network.

7.1. Replacement of filter insert

The equipment is designed and made with the intention to minimize the need of its maintenance. In the regular operation, the replacement interval is once a year - unscrew the filter stopper in the upper part of the chamber, change the filter insert for a new one, screw on the stopper onto the previous place and tighten by hand.

8. TROUBLESHOOTING

W a r n i n g !

All repairs may only be carried out by trained service technicians from competent organizations having an appropriate authorization for such activities.

Before making any interventions into the equipment, it is necessary to reduce the air pressure in the air chamber to zero and to disconnect the equipment from the electric network.

THE DRYING UNIT DOES NOT WORK:

- low operating pressure
 - reduce the air offtake, check the source efficiency, check possible leaks in the distribution
- inoperative solenoid valve
 - repair or change the valve
- clogged nozzle of the regenerating air
 - clean or replace the nozzle, use correct size of the nozzle (see Product maintenance)
- inoperative cooler ventilator
 - replace ventilator
- leak of white liquid through the solenoid valve
 - dismantle the chamber, replace drying substance, bottom filter and sealing and inspect tightness, apply soapy water onto O-rings of the nuts

DRYING UNIT IS NOISY OR SOUNDS ANNOYING :

- defective solenoid valve
 - replace the valve
- damaged damping substance in the condensate vessel
 - replace damping substance or the vessel
- damaged pressure hose
 - replace pressure hose

After troubleshooting and reassembly of the drier, it is necessary to carry out the drier regeneration by starting the compressor with a small air offtake (compressor working mode 20-30%) up to the automatic stop of regeneration.


W a r n i n g !

The manufacturer reserves the right to make modifications to this apparatus without substantially influencing the features thereof.


9. CONNECTION DIAGRAM

1. cooler with ventilator
2. filter stopper
3. humidity sensor chamber
4. return valve
5. solenoid valve of regeneration
6. air chamber
7. adsorption chamber
8. solenoid valve
9. condensate vessel
10. compressor aggregate

10. ELECTRIC WIRING DIAGRAMS

1/N/PE ~ 230 V 50..60 Hz
 1/N/PE ~ 110 V 60 Hz
 ELECTRIC OBJECT CLASS 1
 TYPE B


Legend:

- M Electric motor of compressor 1f., 550W
- EV1 Cooling ventilator of compressor
- EV2 Cooling ventilator of drier
- SP Pressure-operated switch, MDR 2
- ST Thermoswitch, 2450RH
- YV1 Solenoid valve of drier
- YV2 Solenoid valve of regeneration
- C_B Condenser
- FA Safety limit switch 2 pole, D
- SRH Humidity sensor

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Certificate:

Product is manufactured according to production documentation and approved technical conditions. The product is complete and meets all prescribed conditions.

Guarantee:

Manufacturer provides a guarantee of the product in accordance with the paragraphs No.429, 430 and 431 of the Commercial Act No.513/1991, in duration of 24 months from the selling day, when the terms of guarantee document are followed.

- in case of the claim for guarantee repairment a correct fulfilled guarantee document and certificate of quality and completeness of product must be submitted together with the product itself.
- at the guarantee period all failures appeared on the product due to production fault will be removed without any payment.
- the guarantee does not apply to the product damaged during transportation or during incorrect storage.
- the guarantee does not apply to the failures due to wrong operation.
- the guarantee does not apply to the failures due to use of product for any other purposes than are described in the directions for use or as has been understood with the manufacturer.
- the guarantee does not apply to the product when an intervention or arbitrary correction has been made.
- the guarantee does not apply to the product with any completeness which should be discovered during the sale.

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